

Anesthesia-Rx



Key Benefits

- Providers choose Anesthesia-Rx
- Complements existing OR workflow
- Fast and reliable access to medication
- Improves synergy between Pharmacy and OR
- Increases safety and compliance

Providers Choose Anesthesia-Rx for a Fast and Intuitive OR Workflow

Anesthesia-Rx™ is a welcome addition to the operating room. Designed in collaboration with anesthesia providers to support fast and intuitive access to medication, the cart improves communication with pharmacy, enhancing security and regulatory compliance. Anesthesia providers will experience a cart that complements their existing workflow and provides fast and reliable access to medication.

Reliable Access to Medications

Anesthesia-Rx offers providers multiple options to quickly access medications in the event of an emergency.

The **Manual Lock** allows providers to choose when to lock the cart. This feature has been designed so that the cart will not lock after being idle during a procedure, so providers can get the medications they need.

The **Manual Override Button** gives instant access to every pocket in the cart. One push of this button unlocks all of the pockets enabling providers to access what they need during any emergency.

A traditional **Key Lock** is standard on the front of every cart, offering a backup option for accessing medications in the event that all other emergency access features are unavailable.

The **Uninterrupted Power Supply (UPS)** maintains access to medication if power to the cart is lost, with no interruption in delivery of quality patient care, even if the cart is accidentally unplugged or there is a power outage.

Designed to Complement the OR Workflow While Meeting the Needs of Pharmacy

Pharmacy and nursing have been realizing the benefits of automated medication dispensing for more than two decades, but the OR has lagged behind in adoption. Gaining acceptance from anesthesia

providers can prove to be a challenge when implementing automation in the OR. McKesson developed Anesthesia-Rx in conjunction with Anesthesiologists and CRNAs resulting in a cart that providers want to bring into the OR.

RapiDispense™ puts medication in the hands of anesthesia providers, FAST, by simply taking medication out of a non-locked pocket and scanning on dispense.

Multi Touch Dispensing reduces the time spent interacting with the cart allowing providers to focus on the patient. All the medications that need to be dispensed from a drawer can be selected and removed at one time. No need to select and remove medications one at a time, eliminating slow dispensing procedure.

The Virtual Pocket™ Display allows providers to easily find medications on the screen and in the cart. There's no need to search through a list of medications to find what is needed.

MyStock™ provides tracking and secure storage for prepared emergency medication that may be used for a patient at a later time. Medications can be transferred from MyStock and charged to a different patient or another user's stock.

Next Case Preparation supports anesthesia providers preparing for their next case. Toggle between the current case and the next case so that medication can be prepared, tracked, stored and organized securely.

The **Integrated Color Label Printer** reduces the time anesthesia providers spend preparing medication, and helps improve regulatory compliance and patient safety. Every time a medication is dispensed, a color label is automatically printed with all of the necessary information to meet regulatory guidelines. The label is also color coded for easy identification during procedures in order to reduce medication identity confusion.

Expiration Date Tracking provides the ability to remove expired medications from the cart.

Scanning on Restock helps improve patient safety and satisfaction among anesthesia providers. Keeping the right medication in the right pocket reduces the chance for errors and improves safety and organization.

Real-Time Reporting enhances communication between pharmacy and the OR. As medication is stocked, dispensed, returned and wasted, pharmacy is able to view reports to efficiently track inventory.



Anesthesia-Rx is a Class 1 Medical Device and adheres to strict quality standards as defined by the FDA.

MCKESSON

McKesson Automation Solutions

500 Cranberry Woods Drive
Cranberry Township, PA 16066
(724) 741-8000
www.mckesson.com